

PRODUCT INFORMATION

DyLight™ 488-Conjugated Goat Anti-Rabbit IgG Secondary Antibody Item No. 10011231

Overview

Contents:	This vial contains 100 µl of DyLight™ 488-Conjugated Goat anti-Rabbit IgG Secondary Antibody in 50% glycerol. DyLight™ 488 is a product of Thermo Fisher Scientific Inc.
Storage:	-20°C
Stability:	≥2 years

Laboratory Procedures

To prepare a working solution of DyLight™ 488-Conjugated Goat Anti-Rabbit IgG Secondary Antibody, dilute 1:100 in TBS, pH 7.4, containing 0.1% triton X-100. Add 75 to 100 µl of the working solution to each well of your plate and incubate at room temperature for one hour. Longer incubation times may be needed to obtain optimal staining.

Precaution

This product may not perform as described if used with other assay reagents or protocols.

WARNING

THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the [complete](#) Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 01/20/2026

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA

PHONE: [800] 364-9897
[734] 971-3335

FAX: [734] 971-3640

CUSTSERV@CAYMANCHEM.COM
WWW.CAYMANCHEM.COM