


PRODUCT INFORMATION



para-bromo Methylaminorex

Item No. 32835

Formal Name: 5-(4-bromophenyl)-4-methyl-4,5-dihydrooxazol-2-amine
Synonyms: 4'-bromo MAR, *p*-bromo MAR, *para*-bromo MAR, 4'-bromo MAX, *para*-bromo MAX, 4'-bromo Methylaminorex, *p*-bromo Methylaminorex
MF: C₁₀H₁₁BrN₂O
FW: 255.1
Purity: ≥98% (mixture of diastereomers)
Supplied as: A solid
Storage: -20°C
Stability: ≥4 years



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Description

para-bromo Methylaminorex (Item No. 32835) is an analytical reference standard that is structurally similar to known stimulants. This product is intended for research and forensic applications.

WARNING

THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the [complete](#) Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 10/09/2024

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA

PHONE: [800] 364-9897
[734] 971-3335

FAX: [734] 971-3640

CUSTSERV@CAYMANCHEM.COM
WWW.CAYMANCHEM.COM