

PRODUCT INFORMATION

iFluor® 647 Tetrazine

Item No. 44071

FW: 1,272.4
Ex./Em. Max: 656/670 nm
Supplied as: A solid
Storage: -20°C
Stability: ≥4 years
Special Conditions: Light sensitive

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

iFluor® 647 tetrazine is supplied as a solid. A stock solution may be made by dissolving the iFluor® 647 tetrazine in the solvent of choice, which should be purged with an inert gas. iFluor® 647 tetrazine is soluble in DMSO.

Description

iFluor® 647 tetrazine is a fluorophore functionalized with tetrazine, which reacts with *trans*-cyclooctene (TCO) to form a stable conjugate. It displays excitation/emission maxima of 656/670 nm, respectively. iFluor® 647 tetrazine can be used for labeling TCO-modified biological molecules for fluorescence imaging and other fluorescence-based biological applications.

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY
Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 12/02/2025

CAYMAN CHEMICAL
1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA
PHONE: [800] 364-9897
[734] 971-3335
FAX: [734] 971-3640
CUSTSERV@CAYMANCHEM.COM
WWW.CAYMANCHEM.COM