PRODUCT INFORMATION

SM Color Detector

Item No. 10010073

Overview

Contents: This vial contains a lyophilized powder of DAOS (N-ethyl-N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline) and 4-aminoantipyrine.
Storage: 4°C
Stability: ≥2 years

Laboratory Procedures

Prior to use in the Sphingomyelin Colorimetric Assay Kit (Item No. 10009928), reconstitute the Color Detector with 3 ml of diluted SM Buffer (10X) (Item No. 10010072). The reconstituted Color Detector is stable for eight hours at room temperature.

Precaution

This product may not perform as described if used with other assay reagents or protocols.

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent *via* email to your institution.

WARRANTY AND LIMITATION OF REMEDY
Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 01/19/2026

CAYMAN CHEMICAL
1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA
PHONE: [800] 364-9897
[734] 971-3335
FAX: [734] 971-3640
CUSTSERV@CAYMANCHEM.COM
WWW.CAYMANCHEM.COM