
PRODUCT INFORMATION

2’,3’-Dideoxyuridine-5’-O-(1-thiotriphosphate) (sodium salt)

Item No. 38464

Synonyms: ddUTP-α-S, α-S-ddUTP
MF: $C_9H_{11}N_2O_{12}P_3S \bullet 4Na$
FW: 556.1
Purity: ≥95% (a mixture of Rp and Sp isomers)
Supplied as: A solution in water
Storage: -80°C
Stability: ≥4 years



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Description

2’,3’-Dideoxyuridine-5’-O-(1-thiotriphosphate) (ddUTP-α-S) is a sulfur-containing nucleoside triphosphate derivative. ddUTP-α-S is a mixture of Rp and Sp isomers.

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent *via* email to your institution.

WARRANTY AND LIMITATION OF REMEDY
Buyer agrees to purchase the material subject to Cayman’s Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 04/17/2026

CAYMAN CHEMICAL
1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA
PHONE: [800] 364-9897
[734] 971-3335
FAX: [734] 971-3640
CUSTSERV@CAYMANCHEM.COM
WWW.**CAYMANCHEM**.COM