PRODUCT INFORMATION

Immunoassay Buffer B Concentrate (10X)

Item No. 400054

Overview

Contents: This vial contains 10 ml of concentrated Assay Buffer (1 M phosphate, containing 11% BSA, 4 M sodium chloride, 10 mM EDTA, and 0.1% sodium azide).

Stability: ≥3 years at room temperature

Laboratory Procedures

This plate buffer been formulated and tested to work with Cayman's ELISA Kits. Detailed instructions for its use are contained in the respective Cayman ELISA Kit booklet.

Precaution

This product may not perform as described if used with other assay reagents or protocols.

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent *via* email to your institution.

WARRANTY AND LIMITATION OF REMEDY
Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 01/08/2026

CAYMAN CHEMICAL
1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA
PHONE: [800] 364-9897
[734] 971-3335
FAX: [734] 971-3640
CUSTSERV@CAYMANCHEM.COM
WWW.CAYMANCHEM.COM