PRODUCT INFORMATION

NBD Standard

Item No. 401031

Overview

Contents: This vial contains 100 µl of 2-NBDG in ethanol.
Storage: -20°C
Stability: ≥1 year

Laboratory Procedures

NBD Standard has been tested and formulated to work exclusively with Cayman's Myeloperoxidase Chlorination and Peroxidation Inhibitor Screening Assay Kit (Item No. 702930). Please visit www.caymanchem.com/catalog/702930 for the kit protocol, procedures, and product handling.

Precaution

This product may not perform as described if used with other assay reagents or protocols.

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent *via* email to your institution.

WARRANTY AND LIMITATION OF REMEDY
Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 05/13/2025

CAYMAN CHEMICAL
1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA
PHONE: [800] 364-9897
[734] 971-3335
FAX: [734] 971-3640
CUSTSERV@CAYMANCHEM.COM
WWW.**CAYMANCHEM**.COM