PRODUCT INFORMATION

DNA/RNA Oxidative Damage (Clone 7E6.9) ELISA Standard

Item No. 401134

Overview

Contents: This vial contains a 2 µg/ml solution of 8-hydroxy-2'-deoxyguanosine in ethanol.
Storage: -20°C
Stability: ≥2 years

Laboratory Procedures

DNA/RNA Oxidative Damage (Clone 7E6.9) ELISA Standard has been tested and formulated to work exclusively with Cayman's DNA/RNA Oxidative Damage (Clone 7E6.9) ELISA Kit (Item No. 501130). Please visit www.caymanchem.com/catalog/501130 for the kit protocol, procedures, and product handling.

Cautions

This product may not perform as described if used with other assay reagents or protocols.

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent *via* email to your institution.

WARRANTY AND LIMITATION OF REMEDY
Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 03/17/2025

CAYMAN CHEMICAL
1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA
PHONE: [800] 364-9897
[734] 971-3335
FAX: [734] 971-3640
CUSTSERV@CAYMANCHEM.COM
WWW.**CAYMANCHEM**.COM