
PRODUCT INFORMATION

Oligopeptide-68 (acetate)

Item No. 45960



Formal Name:	glycyl-L-arginylglycyl-L-aspartyl-L-tyrosyl-L-isoleucyl-L-tryptophyl-L-seryl-L-leucyl-L-aspartyl-L-threonyl-L-glutamine, acetate
Peptide Sequence:	GRGDYIWSLDTQ-OH
MF:	$C_{62}H_{91}N_{17}O_{21} \bullet XC_2H_4O_2$
FW:	1,410.5
Purity:	≥95%
Supplied as:	A solid
Storage:	-20°C
Stability:	≥4 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Oligopeptide-68 (acetate) is supplied as a solid. A stock solution may be made by dissolving the oligopeptide-68 (acetate) in the solvent of choice, which should be purged with an inert gas. Oligopeptide-68 (acetate) is sparingly soluble (1-10 mg/ml) in acetonitrile:water (1:1).

Description

Oligopeptide-68 is a 12-amino acid peptide, where the first amino acid and last eight amino acids correspond to residues 328-335 of human TGF-β1.

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent *via* email to your institution.

WARRANTY AND LIMITATION OF REMEDY
Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 08/12/2026

CAYMAN CHEMICAL
1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA
PHONE: [800] 364-9897
[734] 971-3335
FAX: [734] 971-3640
CUSTSERV@CAYMANCHEM.COM
WWW.**CAYMANCHEM**.COM