
# PRODUCT INFORMATION

## PROTAC MALT1 Degrader-1

*Item No. 46386*

**CAS Registry No.:** 2975664-37-8
**Formal Name:** N-[5-chloro-6-[5-[[trans-4-[[4-[2-(2,6-dioxo-3-piperidinyl)-2,3 dihydro-1-oxo-1H-isoindol-5-yl]-1-piperazinyl]methyl]cyclohexyl]methyl]-1,2,4-oxadiazol-3-yl]-3-pyridinyl]-N′-[7-(1-methylethyl) pyrazolo[1,5-a]pyrimidin-6-yl]-urea
**Synonyms:** PROTAC Mucosa-associated Lymphoid Tissue Lymphoma Translocation Protein 1 Degrader-1, Proteolysis-targeting Chimera Mucosa-associated Lymphoid Tissue Lymphoma Translocation Protein 1 Degrader-1, Proteolysis-targeting Chimera MALT1 Degrader-1
**MF:** $C_{42}H_{47}ClN_{12}O_5$
**FW:** 835.4
**Purity:** ≥98%
**Supplied as:** A crystalline solid
**Storage:** -20°C
**Stability:** ≥4 years



*Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.*

### Laboratory Procedures

PROTAC MALT1 degrader-1 is supplied as a solid. A stock solution may be made by dissolving the PROTAC MALT1 degrader-1 in the solvent of choice, which should be purged with an inert gas. PROTAC MALT1 degrader-1 is sparingly soluble (1-10 mg/ml) in DMSO and slightly soluble (0.1-1 mg/ml) in a 1:1 solution of acetonitrile:water.

### Description

PROTAC MALT1 degrader-1 is a heterobifunctional molecule that contains a derivative of the mucosa-associated lymphoid tissue lymphoma translocation protein 1 (MALT1) inhibitor MLT-985 conjugated to the cereblon (CRBN) E3 ligase ligand lenalidomide (Item No. 14643) *via* a linker.[1]

### Reference


1. Zhang, R.Y., Wang, Z.X., Zhang, M.Y., *et al*. MALT1 inhibitors and degraders: Strategies for NF-κB-driven malignancies. *J. Med. Chem.* **68(5)**, 5075-5096 (2025).



**WARNING**
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

**SAFETY DATA**
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent *via* email to your institution.

**WARRANTY AND LIMITATION OF REMEDY**
Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 09/09/2026



**CAYMAN CHEMICAL**
1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA
**PHONE:** [800] 364-9897
[734] 971-3335
**FAX:** [734] 971-3640
CUSTSERV@CAYMANCHEM.COM
WWW.**CAYMANCHEM**.COM