PRODUCT INFORMATION

Glycolysis Assay Cofactor

Item No. 600454

Overview

Contents: This vial contains 250 µl of a 100X cofactor solution.
Storage: -20°C
Stability: ≥2 years

Laboratory Procedures

Glycolysis assay cofactor has been tested and formulated to work exclusively with Cayman's Glycolysis Cell-Based Assay Kit and Oxygen Consumption/Glycolysis Dual Assay Kit (Item Nos. 600450 and 601060, respectively). Please visit www.caymanchem.com/catalog/600450 or www.caymanchem.com/catalog/601060 for the kit protocol, procedures, and product handling.

Precaution

This product may not perform as described if used with other assay reagents or protocols.

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent *via* email to your institution.

WARRANTY AND LIMITATION OF REMEDY
Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 11/20/2025

CAYMAN CHEMICAL
1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA
PHONE: [800] 364-9897
[734] 971-3335
FAX: [734] 971-3640
CUSTSERV@CAYMANCHEM.COM
WWW.**CAYMANCHEM**.COM