PRODUCT INFORMATION

Uric Acid Enzyme Mixture

Item No. 700323

Overview

Contents: This vial contains a lyophilized powder of uricase and horseradish peroxidase.
Storage: -20°C
Stability: ≥2 years

Laboratory Procedures

Prior to use in the Uric Acid Assay Kit (Item No. 700320), add 900 µl of 1X Uric Acid Assay Buffer (Item No. 700321) to the vial, vortex, and place the vial on ice. This is enough enzyme mixture to assay 60 wells. The reconstituted mixture is stable for four hours on ice. Please visit www.caymanchem.com/catalog/kit number 700320 for the kit protocol, procedures, and product handling.

Precaution

This product may not perform as described if used with other assay reagents or protocols.

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent *via* email to your institution.

WARRANTY AND LIMITATION OF REMEDY
Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 01/26/2026

CAYMAN CHEMICAL
1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA
PHONE: [800] 364-9897
[734] 971-3335
FAX: [734] 971-3640
CUSTSERV@CAYMANCHEM.COM
WWW.**CAYMANCHEM**.COM